

25 ml Conical Tube

(#50025/50225/54025/50125)

1. Description

SPL 25 ml Conical Tube is optimized for storage of samples with a volume of 25 ml or less. The shorter tube compared to the 50 ml tube is easy to use and minimizes sample contamination. It can be used with a 50 ml centrifuge rotor using the dedicated adapter.

2. Features

- ❖ Tight sealing screw cap
- ❖ Imprinted external graduations with marking area
- ❖ Packed in zipper bags
- ❖ 25 ml Tube Adapter for 50 ml rotor (Cat. No. 52025)
- ❖ Amber Type (Cat. No. 54025)
- ❖ Provided in racks (Cat. No. 50125)

3. Functional Test

Non-Cytotoxic	:	Passed
Non-Pyrogenic	:	Passed
DNase/RNase Free	:	Passed
Human DNA free	:	Passed

4. Product Range

In Sleeves

Cat. No.	Material (Tube / Cap)	Color	Bottom Type	External Dimensions d x h (mm)	Total Vol. (ml)	RCF Rating	Sterile	Packaging
50025	PP / HDPE	Clear	Conical	30.00 x 77.20	25.00	18,000 x g	+	25 / 200
50225	PP / HDPE	Clear	Skirted	30.00 x 77.20	25.00	-	+	25 / 200
54025	PP / HDPE	Amber	Conical	30.00 x 77.20	25.00	18,000 x g	+	25 / 200

In Racks

Cat. No.	Material (Tube / Cap)	Color	Bottom Type	External Dimensions d x h (mm)	Total Vol. (ml)	RCF Rating	Sterile	Packaging
50125	PP / HDPE	Clear	Conical	30.00 x 77.20	25.00	18,000 x g	+	25 / 200

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs



SPL-TDS-CRY




7. Specifications & Technical Drawing

(1) Conical type (#50025/54025)
- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 34.04 mm Height: 13.00 mm	High-density polyethylene (HDPE)	<ul style="list-style-type: none"> • Pieces/Sleeve: 25 • Pieces/Case: 200
Body	Diameter: 30.00 mm Height: 77.20 mm	Polypropylene (PP)	


- Technical Drawing



(2) Skirted type (#50225)
- Specifications

	nDimensions	Material (Resin)	Packaging
Cap	Diameter: 34.04 mm Height: 13.00 mm	High-density polyethylene (HDPE)	<ul style="list-style-type: none"> • Pieces/Sleeve: 25 • Pieces/Case: 200
Body	Diameter: 30.00 mm Height: 78.20 mm	Polypropylene (PP)	

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
To place an order, contact your local distributor or
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-CRY

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 4, Nov 7th, 2019