

## 5 ml Snap Tube & Adapter

(#50005/50105/55005/55105/52005)

### 1. Description

5 ml Snap tubes conical type are widely used in laboratories for multiple purposes, from storage of samples to high speed centrifugation. The skirt type is self-standing, so it is convenient to handle. Snap tubes can assist users during experiments by enabling easy one-hand opening and closing of the tubes by providing snap type cap.

### 2. Features

- ❖ Snap cap with tight sealing
- ❖ Imprinted external graduations with marking area
- ❖ Sterile (Cat. No.50105, 55105)
- ❖ Skirted, self-standing type (Cat. No. 55005, 55105)
- ❖ 5 ml Snap Tube Adapter for 15 ml rotor (Cat.No.52005)

### 3. Functional Test

Sterility	:	Passed
Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase Free	:	Passed
DNA free	:	Passed

### 4. Product Range

Cat.No.	Type	Material	External Dimension (mm)	Total Volume (ml)	RCF Rating	Sterile	Packaging
50005	Conical	PP	20.00 x 62.00	5.00	25,000 x g	-	100 in sleeve / 200
50105	Conical	PP	20.00 x 62.00	5.00	25,000 x g	+	100 in sleeve / 200
55005	Skirted	PP	20.00 x 63.30	5.00	-	-	100 in sleeve / 200
55105	Skirted	PP	20.00 x 63.30	5.00	-	+	100 in sleeve / 200

Cat.No.	Material	External Dimension (mm)	RCF Rating	Packaging
52005	Acetal	16.60 x 78.24	25,000 x g	10

### 5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.


## 6. Specifications & Technical Drawing

(1) 5 ml Snap Tube (#50005/#50105)

- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 22.08 mm Height: 7.00 mm	Polypropylene(PP)	<ul style="list-style-type: none"> <li>▪ Pieces/sleeve:100</li> <li>▪ Pieces/case: 200</li> </ul>
Body	Diameter: 15.80 mm Height: 62.00 mm Volume: 5.00 ml	Polypropylene(PP)	

- Technical Drawing




(2) 5 ml Snap Tube, skirted (#55005/#55105)

- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 22.08 mm Height: 7.00 mm	Polypropylene(PP)	<ul style="list-style-type: none"> <li>▪ Pieces/sleeve:100</li> <li>▪ Pieces/case: 200</li> </ul>
Body	Diameter: 15.19 mm Height: 63.30 mm Volume: 5.00 ml	Polypropylene(PP)	

- Technical Drawing



SPL-TDS-5ST


For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.  
© SPL Life Sciences Co., Ltd., Version 4, JAN 11<sup>th</sup>, 2021

## (3) 5 ml Snap Tube Adapter (#52005)

## - Specifications

	Dimensions	Material (Resin)	Packaging
Adapter	Diameter: 16.60 mm Height: 78.24 mm	Acetal	▪ Pieces/case: 10

## - Technical Drawing



[www.spllifesciences.com](http://www.spllifesciences.com)

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [spl@ispl.co.kr](mailto:spl@ispl.co.kr)

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [business@ispl.co.kr](mailto:business@ispl.co.kr)