

Amber Conical Tube (#54015, 54050)

1. Description

Amber Conical Tubes are suitable for handling and storing light-sensitive samples. Materials used for amber conical tubes significantly decrease light transmission over the visible wavelength spectrum compared to clear polypropylene tubes.

Also, Amber tubes are feasible to centrifugation, applicable to wide range of experiments.

2. Features

- ❖ Light protection for light-sensitive samples
- ❖ Screw cap with tight sealing
- ❖ External graduations with marking area
- ❖ Packed in zipper bags

3. Functional Test

Sterility	:	Passed
Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase Free	:	Passed
DNA free	:	Passed

4. Product Range

Cat.No.	Bottom Type	Material (Body/Cap)	External Dimension d x h (mm)	Total Volume (ml)	RCF Rating	Sterile	Packaging
54015	Conical	PP / HDPE	17.00 x 120.70	15.00	13000 x g	+	50 in Sleeve / 500
54115	Conical	PP / HDPE	17.00 x 120.70	15.00	13000 x g	+	50 in Sleeve / 200
54050	Conical	PP / HDPE	30.00 x 116.70	50.00	14000 x g	+	25 in Sleeve / 500
54150	Conical	PP / HDPE	30.00 x 116.70	50.00	14000 x g	+	25 in Sleeve / 200

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs



SPL-TDS-AMCTU

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 3, July 19th, 2016


7. Specifications & Technical Drawing

(1) 15 ml Conical Tube (# 54015, 54115)

- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 21.30 mm Height: 11.00 mm	High-density polyethylene (HDPE)	#54015 ▪ Pieces/sleeve: 50 ▪ Pieces/case: 500
Body	Diameter: 16.90 mm Height: 119.35 mm Volume: 15.00 ml	Polypropylene(PP)	#50115 ▪ Pieces/sleeve: 50 ▪ Pieces/case: 200

- Technical Drawing




(2) 50 ml Conical Tube (# 54050, 54150)

- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 34.00 mm Height: 13.00 mm	High-density polyethylene (HDPE)	#54050 ▪ Pieces/sleeve: 25 ▪ Pieces/case: 500
Body	Diameter: 29.30 mm Height: 114.10 mm Volume: 50.00 ml	Polypropylene(PP)	#54150 ▪ Pieces/sleeve: 25 ▪ Pieces/case: 200

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-AMCTU

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 3, July 19th, 2016